

GlyT2 Rabbit mAb

Catalog No: #52316



Package Size: #52316-1 50ul #52316-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	GlyT2 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S08-1A3
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human Glyt2
Conjugates	Unconjugated
Modification	Unmodification
Other Names	NET1; GLYT2; HKPX3; GLYT-2
Accession No.	Swiss-Prot:Q9Y345GeneID:9152
Uniprot	Q9Y345
GeneID	9152
Calculated MW	Calculated MW: 87 kDa; Observed MW: 87 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000

Images



Western blot detection of Glyt2 in Jurkat,C6 cell lysates using Glyt2 Rabbit mAb(1:1000 diluted).Predicted band size:87kDa.Observed band size:87kDa.

Background

Swiss-Prot Acc.Q9Y345.Sodium- and chloride-dependent glycine transporter 2

Note: This product is for in vitro research use only