

LGALS14 Rabbit mAb

Catalog No: #52396



Package Size: #52396-1 50ul #52396-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	LGALS14 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S09-6C1
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse
Immunogen Description	A synthetic peptide of human LGALS14
Conjugates	Unconjugated
Modification	Unmodification
Other Names	Charcot-Leyden crystal protein 2;CLC2;Galectin-14;Gal-14
Accession No.	Swiss-Prot:Q8TCE9GeneID:56891
Uniprot	Q8TCE9
GeneID	56891
Calculated MW	Calculated MW: 16 kDa; Observed MW: 16 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000;

Images



Western blot detection of LGALS14 in C6 cell lysates using LGALS14 Rabbit mAb(1:1000 diluted).Predicted band size:16kDa.Observed band size:16kDa.

Background

Swiss-Prot Acc.Q8TCE9.Placental protein 13-like

Note: This product is for in vitro research use only