

Prealbumin Rabbit mAb

Catalog No: #52521



Package Size: #52521-1 50ul #52521-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Prealbumin Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S07-4A2
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human, Rat
Immunogen Description	A synthetic peptide of human Prealbumin
Conjugates	Unconjugated
Modification	Unmodification
Other Names	CTS; ATTR; CTS1; PALB; TBPA; HEL111; HsT2651
Accession No.	Swiss-Prot:P02766GenelD:7276
Uniprot	P02766
GenelD	7276
Calculated MW	Calculated MW: 16 kDa; Observed MW: 16 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000;

Images



Western blot detection of Prealbumin in Rat Brain lysates using Prealbumin Rabbit mAb(1:1000 diluted).Predicted band size:16kDa.Observed band size:16kDa.

Background

Swiss-Prot Acc.P02766.Thyroid hormone-binding protein. Probably transports thyroxine from the bloodstream to the brain.

Note: This product is for in vitro research use only