

SAMD9 Rabbit mAb

Catalog No: #52566



Package Size: #52566-1 50ul #52566-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	SAMD9 Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S09-5B4
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human
Immunogen Description	A synthetic peptide of human SAMD9
Conjugates	Unconjugated
Modification	Unmodification
Other Names	SAM domain-containing protein 9
Accession No.	Swiss-Prot:Q5K651GeneID:54809
Uniprot	Q5K651
GeneID	54809
Calculated MW	Calculated MW: 184 kDa; Observed MW: 184 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000;

Images



Western blot detection of SAMD9 in A549 cell lysates using SAMD9 Rabbit mAb(1:1000 diluted).Predicted band size:184kDa.Observed band size:184kDa.

Background

Swiss-Prot Acc.Q5K651.Sterile alpha motif domain-containing protein 9

Note: This product is for in vitro research use only