

Translin Associated Protein X Rabbit mAb

Catalog No: #52629



Package Size: #52629-1 50ul #52629-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


Description

Product Name	Translin Associated Protein X Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S02-3D3
Isotype	Rabbit IgG
Purification	Affinity Purified
Applications	WB
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human TRAX
Conjugates	Unconjugated
Modification	Unmodification
Other Names	TRAX
Accession No.	Swiss-Prot:Q99598GeneID:7257
Uniprot	Q99598
GeneID	7257
Calculated MW	Calculated MW: 33 kDa; Observed MW: 33 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

Application Details

WB: 1/1000-1/5000

Images



Western blot detection of TRAX in K562,Rat Brain,C6,Hela cell lysates using TRAX Rabbit mAb(1:1000 diluted).Predicted band size:33kDa.Observed band size:33kDa.

## Background

---

Swiss-Prot Acc.Q99598.Acts in combination with TSN as an endonuclease involved in the activation of the RNA-induced silencing complex (RISC).  
Possible role in spermatogenesis.

---

Note: This product is for in vitro research use only