

ADK Rabbit mAb

Catalog No: #52711



Package Size: #52711-1 50ul #52711-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ADK Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal antibody
Clone No.	S05-0G0
Isotype	IgG
Purification	Affinity Purified
Applications	WB IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human ADK
Conjugates	Unconjugated
Modification	Unmodification
Other Names	AK
Accession No.	Swiss-Prot:P55263GenelD:132
Uniprot	P55263
GenelD	132
Calculated MW	Calculated MW:41 kDa,Observed MW:41 kDa
Concentration	0.3 mg/ml
Formulation	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.


Application Details

WB: 1/1000 IHC: 1/100

Images



Western blot detection of ADK in Hela,A549,HL-60,U251,U87-MG cell lysates using ADK Rabbit mAb(1:1000 diluted).Predicted band size:41kDa.Observed band size:41kDa.



Immunohistochemistry of ADK in paraffin-embedded Human tonsil using ADK Rabbit mAb at dilution 1/50

Background

ATP dependent phosphorylation of adenosine and other related nucleoside analogs to monophosphate derivatives. Serves as a potential regulator of concentrations of extracellular adenosine and intracellular adenine nucleotides.

Note: This product is for in vitro research use only