

C17orf74 Rabbit Polyclonal Antibody

Catalog No: #53195



Package Size: #53195-1 50ul #53195-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	C17orf74 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthetic peptide of human C17orf74 (NP_783861.3).
Other Names	C17orf74
Accession No.	Swiss Prot:Q0P670Gene ID:201243
Uniprot	Q0P670
GeneID	201243
Calculated MW	56kDa
SDS-PAGE MW	65kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200

Images



Western blot analysis of extracts of various cell lines, using C17orf74 antibody.



Immunofluorescence analysis of NIH3T3 cells using C17orf74 Rabbit pAb.

Background

—

Note: This product is for in vitro research use only