

C5orf52 Rabbit Polyclonal Antibody

Catalog No: #53295



Package Size: #53295-1 50ul #53295-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	C5orf52 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of mouse C5orf52 (NP_080426.1).
Other Names	C5orf52;4921536K21Rik;AV043724
Accession No.	Swiss Prot:Q9CR34Gene ID:67430
Uniprot	Q9CR34
GeneID	67430
Calculated MW	10kDa/19kDa
SDS-PAGE MW	18kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:2000 - 1:5000IHC 1:50 - 1:200


Images




Western blot analysis of extracts of various cell lines, using C5orf52 antibody.



Immunohistochemistry of paraffin-embedded rat brain using C5orf52 Rabbit pAb.



Immunohistochemistry of paraffin-embedded human placenta using C5orf52 Rabbit pAb.



Immunohistochemistry of paraffin-embedded mouse brain using C5orf52 Rabbit pAb.

Background

Note: This product is for in vitro research use only