

PAXIP1 Rabbit Polyclonal Antibody

Catalog No: #54417



Package Size: #54417-1 50ul #54417-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PAXIP1 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human PAXIP1.
Other Names	PAXIP1;CAGF28;CAGF29;PACIP1;PAXIP1L;PTIP;TNRC2
Accession No.	Uniprot:Q6ZW49GeneID:22976
Uniprot	Q6ZW49
GeneID	22976
Calculated MW	93kDa/117kDa/121kDa
SDS-PAGE MW	121KDa
Formulation	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB1:500 - 1:1000IHC1:50 - 1:200


Images




Western blot analysis of extracts of various cell lines, using PAXIP1 antibody.




Immunohistochemistry of paraffin-embedded rat ovary using PAXIP1 Rabbit pAb.




Immunohistochemistry of paraffin-embedded mouse brain using PAXIP1 Rabbit pAb.




Immunohistochemistry of paraffin-embedded rat spleen using PAXIP1 Rabbit pAb.



Immunohistochemistry of paraffin-embedded human breast cancer using PAXIP1 Rabbit pAb.



Immunohistochemistry of paraffin-embedded mouse spleen using PAXIP1 Rabbit pAb.



Immunohistochemistry of paraffin-embedded human thyroid cancer using PAXIP1 Rabbit pAb.

Background

This gene is a member of the paired box (PAX) gene family and encodes a nuclear protein with six BRCT (breast cancer carboxy-terminal) domains. This protein plays a critical role in maintaining genome stability, condensation of chromatin and progression through mitosis.

Note: This product is for in vitro research use only