

AP3S2 Rabbit Polyclonal Antibody

Catalog No: #54439



Package Size: #54439-1 50ul #54439-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	AP3S2 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant fusion protein of human AP3S2 (NP_005820.1).
Other Names	AP3S2;AP3S3;sigma3b
Accession No.	Swiss Prot:P59780Gene ID:10239
Uniprot	P59780
GeneID	10239
Calculated MW	14kDa/22kDa/43kDa
SDS-PAGE MW	22kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IF 1:50 - 1:200


Images



Western blot analysis of extracts of various cell lines, using AP3S2 Rabbit pAb.



Immunofluorescence analysis of C6 cells using AP3S2 antibody.



Immunofluorescence analysis of U-2 OS cells using AP3S2 antibody.

Background

—

Note: This product is for in vitro research use only