

ASB3 Rabbit Polyclonal Antibody

Catalog No: #54577



Package Size: #54577-1 50ul #54577-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	ASB3 Rabbit Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Affinity purification
Applications	WB,IHC,IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	Recombinant protein of human ASB3.
Other Names	—
Accession No.	Swiss Prot:Q9Y575Gene ID:51130
Uniprot	Q9Y575
GeneID	51130
Calculated MW	—
SDS-PAGE MW	58kDa
Formulation	Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.


Application Details

WB 1:500 - 1:2000IHC 1:50 - 1:200IF 1:50 - 1:200


Images




Western blot analysis of extracts of various cell lines, using ASB3 antibody.




Immunohistochemistry of paraffin-embedded rat kidney using ASB3 Rabbit pAb.




Immunohistochemistry of paraffin-embedded human colon carcinoma using ASB3 Rabbit pAb.



Immunohistochemistry of paraffin-embedded mouse liver using ASB3 Rabbit pAb.



Immunofluorescence analysis of NIH-3T3 cells using ASB3 Rabbit pAb.



Immunofluorescence analysis of C6 cells using ASB3 Rabbit pAb.

—

Note: This product is for in vitro research use only