

## KGF Antibody FITC Conjugated

Catalog No: #C00415F

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	KGF Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 150-194 194 derived from human KGF Fibroblast growth factor 7
Conjugates	FITC
Target Name	KGF
Other Names	KGF; HBGF-7; Fibroblast growth factor 7; FGF-7; Heparin-binding growth factor 7; Keratinocyte growth factor; FGF7
Accession No.	Swiss-Prot#P21781NCBI Gene ID2252
Uniprot	P21781
GeneID	2252;
Excitation Emission	494nm 518nm
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

## Background

Plays an important role in the regulation of embryonic development, cell proliferation and cell differentiation. Required for normal branching morphogenesis. Growth factor active on keratinocytes. Possible major paracrine effector of normal epithelial cell proliferation.

Note: This product is for in vitro research use only