

CDKN1B Antibody PE Conjugated

Catalog No: #C00421P

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CDKN1B Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 150-198 198 derived from human CDKN1B
Conjugates	PE
Target Name	CDKN1B
Other Names	KIP1; MEN4; CDKN4; MEN1B; P27KIP1; Cyclin-dependent kinase inhibitor 1B; Cyclin-dependent kinase inhibitor p27; CDKN1B
Accession No.	Swiss-Prot#P46527NCBI Gene ID1027
Uniprot	P46527
GeneID	1027;
Excitation Emission	480,565nm 578nm
Cell Localization	Cytoplasm, Nucleus
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Background

Important regulator of cell cycle progression. Involved in G1 arrest. Potent inhibitor of cyclin E- and cyclin A-CDK2 complexes. Forms a complex with cyclin type D-CDK4 complexes and is involved in the assembly, stability, and modulation of CCND1-CDK4 complex activation. Acts either as an inhibitor or an activator of cyclin type D-CDK4 complexes depending on its phosphorylation state and or stoichiometry.

Note: This product is for in vitro research use only