

Lysozyme Antibody FITC Conjugated

Catalog No: #C00458F

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lysozyme Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	Full length native Lysozyme
Conjugates	FITC
Target Name	Lysozyme
Other Names	LZM; Lysozyme C; 1,4-beta-N-acetylmuramidase C; LYZ
Accession No.	Swiss-Prot#P61626NCBI Gene ID4069
Uniprot	P61626
GeneID	4069;
Excitation Emission	494nm 518nm
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Lysozymes have primarily a bacteriolytic function; those in tissues and body fluids are associated with the monocyte-macrophage system and enhance the activity of immunoagents.

Note: This product is for in vitro research use only