

HHV8 ORF50 Antibody HRP Conjugated

Catalog No: #C00473H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HHV8 ORF50 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	herpesvirus8typeP
Crossing Reactivity	HHV8
Immunogen Description	KLH conjugated synthetic peptide aa 500-550 691 derived from human HHV8 ORF50
Conjugates	HRP
Target Name	HHV8 ORF50
Other Names	HHV8 ORF50; HHV 8; Human Herpes Virus 8; Human herpesvirus 8; Kaposi's sarcoma associated herpes virus; Kaposi's sarcoma associated herpesvirus; KSHV.
Excitation Emission	N A
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

HHV8 has been found to be associated with three different diseases observed in AIDS patients; kaposi's sarcoma, primary effusion lymphoma (which is a rare type of non-Hodgkin lymphoma affecting the body cavities) and multicentric Castleman's disease. To date there is much evidence to support a direct role for HHV8 in kaposi's sarcoma.

Note: This product is for in vitro research use only