

## FSH receptor Antibody FITC Conjugated

Catalog No: #C00485F

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	FSH receptor Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 225-270 695 derived from human FSHR
Conjugates	FITC
Target Name	FSH receptor
Other Names	LGR1; ODG1; FSHRO; Follicle-stimulating hormone receptor; FSH-R; Folitropin receptor; FSHR
Accession No.	Swiss-Prot#P23945NCBI Gene ID2492
Uniprot	P23945
GeneID	2492;
Excitation Emission	494nm 518nm
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

## Background

Receptor for follicle-stimulating hormone. The activity of this receptor is mediated by G proteins which activate adenylate cyclase.

Note: This product is for in vitro research use only