

FSH receptor Antibody HRP Conjugated

Catalog No: #C00485H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FSH receptor Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 225-270 695 derived from human FSHR
Conjugates	HRP
Target Name	FSH receptor
Other Names	LGR1; ODG1; FSHRO; Follicle-stimulating hormone receptor; FSH-R; Folitropin receptor; FSHR
Accession No.	Swiss-Prot#P23945NCBI Gene ID2492
Uniprot	P23945
GeneID	2492;
Excitation Emission	N A
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Receptor for follicle-stimulating hormone. The activity of this receptor is mediated by G proteins which activate adenylate cyclase.

Note: This product is for in vitro research use only