

SCRN3 Antibody

Catalog No: #25156




Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	SCRN3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity chromatography purified via peptide column
Applications	ELISA WB ICC
Species Reactivity	Hu Ms
Immunogen Type	Peptide
Immunogen Description	Raised against a 15 amino acid peptide near the center of human SCRIN3.
Target Name	SCRIN3
Other Names	Secernin 3, SES3
Accession No.	Swiss-Prot:Q0VDG4Gene ID:79634
Uniprot	Q0VDG4
GeneID	79634;
Concentration	1mg/ml
Formulation	Supplied in PBS containing 0.02% sodium azide.
Storage	Can be stored at -20°C, stable for one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Images



Western blot analysis of SCRIN3 in 293 cell tissue lysate with SCRIN3 antibody at (A) 1 and (B) 2 ug/mL.



Immunocytochemistry of SCRIN3 in 293 cells with SCRIN3 antibody at 20 ug/mL.

Background

SCRN3 is a member of the secernin family of proteins which includes the homologous SCR1, a protein that was first identified as being involved in the regulation of exocytosis from peritoneal mast cells. Little is known of SCR3, but studies have shown that SCR1 may possess epitopes that function as tumor-associated antigens in gastric cancers and increased SCR1 expression in patients with colorectal cancer correlated with poor prognosis, suggesting that SCR3 may also be involved in protein secretion or the regulation of cell growth.

Note: This product is for in vitro research use only