

## IRF-4 Polyclonal Antibody

Catalog No: #44020



Package Size: #44020-1 50ul #44020-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	IRF-4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB,IHC-p,IF/ICC,ELISA
Species Reactivity	Human,Mouse
Specificity	IRF-4 Polyclonal Antibody detects endogenous levels of IRF-4 protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human IRF4.
Accession No.	Swiss-Prot#: Q15306
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western Blot: 1/500 - 1/2000

Immunohistochemistry: 1/100 - 1/300

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/20000.

Note: This product is for in vitro research use only