

## FAS-L Polyclonal Antibody

Catalog No: #44059



Package Size: #44059-1 50ul #44059-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	FAS-L Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Applications	WB;IHC-p;IF/ICC;ELISA
Species Reactivity	Human,Mouse
Specificity	FAS-L Polyclonal Antibody detects endogenous levels of FAS-L protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human FAS ligand.
Accession No.	Swiss-Prot#: P48023 NCBI Gene ID:356
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western Blot: 1/500 - 1/2000

Immunohistochemistry: 1/100 - 1/300

Immunofluorescence: 1/200 - 1/1000.ELISA: 1/40000

Note: This product is for in vitro research use only