

## CNPase Antibody FITC Conjugated

Catalog No: #C00535F

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CNPase Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 145-195 421 derived from human CNPase
Conjugates	FITC
Target Name	CNPase
Other Names	CNP1; 2',3'-cyclic-nucleotide 3'-phosphodiesterase; CNP; CNPase
Accession No.	Swiss-Prot#P09543NCBI Gene ID1267
Uniprot	P09543
GeneID	1267;
Excitation Emission	494nm 518nm
Cell Localization	Extracellular, Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

## Background

May participate in RNA metabolism in the myelinating cell. CNP is the third most abundant protein in central nervous system myelin.

Note: This product is for in vitro research use only