

JNK1/2/3 Rabbit mAb

Catalog No: #56078



Package Size: #56078-1 50ul #56078-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	JNK1/2/3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF IP
Species Reactivity	Human Mouse Rat Cow Dog Monkey Zebrafish
Specificity	JNK1+JNK2+JNK3 Antibody detects endogenous levels of total JNK1+JNK2+JNK3
Immunogen Description	A synthesized peptide derived from human JNK1+JNK2+JNK3
Other Names	c-Jun N-terminal kinase 3; JNK1; JNK2;JNK3; JNK3 alpha protein kinase; JNK3A; MAP kinase 10; MAP kinase p49 3F12; MAPK 10; MAPK10;
Accession No.	Uniprot:P45983/P45984/P53779
Uniprot	P45983/P45984/P53779
Calculated MW	46,54kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.


Application Details

WB:1:1000~1:2000ICC/IF:1:50~1:200IP:1:50

Images



Western blot analysis of JNK1/2/3 Antibody expression in HeLa cell lysate.



Immunofluorescent analysis of HeLa cells, using JNK1/2/3 Antibody.

Product Description

Serine/threonine-protein kinase involved in various processes such as cell proliferation, differentiation, migration, transformation and programmed cell death. Extracellular stimuli such as proinflammatory cytokines or physical stress stimulate the stress-activated protein kinase/c-Jun N-terminal kinase (SAP/JNK) signaling pathway.

Note: This product is for in vitro research use only