

NCR1 Antibody HRP Conjugated

Catalog No: #C00543H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NCR1 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 170-220 304 derived from human CD335
Conjugates	HRP
Target Name	NCR1
Other Names	LY94; CD335; NKP46; NK-p46; Natural cytotoxicity triggering receptor 1; Lymphocyte antigen 94 homolog; NK cell-activating receptor; Natural killer cell p46-related protein; hNKp46; NCR1
Accession No.	Swiss-Prot#O76036NCBI Gene ID9437
Uniprot	O76036
GeneID	9437;
Excitation Emission	N A
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000

Background

Cytotoxicity-activating receptor that may contribute to the increased efficiency of activated natural killer (NK) cells to mediate tumor cell lysis.

Note: This product is for in vitro research use only