

HLA G Rabbit mAb

Catalog No: #56134

Package Size: #56134-1 50ul #56134-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	HLA G Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	HLA G Antibody detects endogenous levels of total HLA G
Immunogen Description	A synthesized peptide derived from human HLA G
Other Names	HLA G antigen; MHC class I antigen G; HLA-G; HLA60; HLAG;
Accession No.	Uniprot:P17693
Uniprot	P17693
Calculated MW	42kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000

Images



Western blot analysis of HLA G expression in Jurkat cell lysate.

Product Description

Involved in the presentation of foreign antigens to the immune system. Plays a role in maternal tolerance of the fetus by mediating protection from the deleterious effects of natural killer cells, cytotoxic T lymphocytes, macrophages and mononuclear cells.

Note: This product is for in vitro research use only