

TMEM16A Rabbit mAb

Catalog No: #56151



Package Size: #56151-1 50ul #56151-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TMEM16A Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	IHC FC
Species Reactivity	Human
Specificity	TMEM16A Antibody detects endogenous levels of total TMEM16A
Immunogen Description	A synthesized peptide derived from human TMEM16A
Other Names	Anoctamin-1; ANO1; DOG1; ORAOV2; TAOS2; TMEM16A;
Accession No.	Uniprot:Q5XXA6
Uniprot	Q5XXA6
Calculated MW	kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

IHC:1:50~1:2000FC:1:50

Product Description

Acts as a calcium-activated chloride channel. Required for normal tracheal development.

Note: This product is for in vitro research use only