

Gata6 Rabbit mAb

Catalog No: #56244

Package Size: #56244-1 50ul #56244-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Gata6 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF
Species Reactivity	Human
Specificity	Gata6 Antibody detects endogenous levels of total Gata 6
Immunogen Description	A synthesized peptide derived from human Gata 6
Other Names	Gata binding factor 6; GATA-binding factor 6; Gata6; Transcription factor Gata 6;
Accession No.	Uniprot:Q92908
Uniprot	Q92908
Calculated MW	60kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000ICC/IF:1:50~1:200

Images



Western blot analysis of Gata 6 expression in human fetal heart lysate.

Product Description

Thought to be important for regulating terminal differentiation and/or proliferation.

Note: This product is for in vitro research use only