

CD97 Rabbit mAb

Catalog No: #56288

Package Size: #56288-1 50ul #56288-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CD97 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human Mouse
Specificity	CD97 Antibody detects endogenous levels of total CD97
Immunogen Description	A synthesized peptide derived from human CD97
Other Names	ADGRE5; Adhesion G protein-coupled receptor E5; CD97 antigen; Leukocyte antigen CD97; Seven span transmembrane protein; TM7LN1;
Accession No.	Uniprot:P48960
Uniprot	P48960
Calculated MW	92kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:50FC:1:50

Images



Western blot analysis of CD97 expression in U937 cell lysate.

Product Description

Receptor potentially involved in both adhesion and signaling processes early after leukocyte activation. Plays an essential role in leukocyte migration.

Note: This product is for in vitro research use only