

GRK3 Rabbit mAb

Catalog No: #56322

Package Size: #56322-1 50ul #56322-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GRK3 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF
Species Reactivity	Human
Specificity	GRK3 Antibody detects endogenous levels of total GRK3
Immunogen Description	A synthesized peptide derived from human GRK3
Other Names	ADRBK2; BARK2; GRK3;
Accession No.	Uniprot:P35626
Uniprot	P35626
Calculated MW	80kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200

Images



Western blot analysis of GRK3 expression in Ramos cell lysate.

Product Description

Specifically phosphorylates the agonist-occupied form of the beta-adrenergic and closely related receptors.

Note: This product is for in vitro research use only