

## Netrin 1 Rabbit mAb

Catalog No: #56336

Package Size: #56336-1 50ul #56336-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	Netrin 1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF
Species Reactivity	Human Mouse Rat
Specificity	Netrin 1 Antibody detects endogenous levels of total Netrin 1
Immunogen Description	A synthesized peptide derived from human Netrin 1
Other Names	Netrin1; NTN1; NTN1L; Unc6;
Accession No.	Uniprot:O95631
Uniprot	O95631
Calculated MW	75kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000ICC/IF:1:50~1:200

## Images



Western blot analysis of Netrin 1 expression in Human fetal brain lysate.

## Product Description

Netrins control guidance of CNS commissural axons and peripheral motor axons. Its association with either DCC or some UNC5 receptors will lead to axon attraction or repulsion, respectively.

---

Note: This product is for in vitro research use only