

TBX21 Rabbit mAb

Catalog No: #56342

Package Size: #56342-1 50ul #56342-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TBX21 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	IHC
Species Reactivity	Human
Specificity	T-bet / Tbx21 Antibody detects endogenous levels of total T-bet / Tbx21
Immunogen Description	A synthesized peptide derived from human T-bet / Tbx21
Other Names	T bet; T box 21; T PET; TBET; TBLYM; Tbx21; TPET;
Accession No.	Uniprot:Q9UL17
Uniprot	Q9UL17
Calculated MW	kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Application Details

IHC:1:50~1:200

Images



Immunohistochemical analysis of paraffin-embedded human spleen, using T-bet/Tbx21 Antibody.

Product Description

Transcription factor that controls the expression of the TH1 cytokine, interferon-gamma. Initiates TH1 lineage development from naive TH precursor cells both by activating TH1 genetic programs and by repressing the opposing TH2 programs.

Note: This product is for in vitro research use only