

Alpha s2 Casein Antibody HRP Conjugated

Catalog No: #C00550H

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Alpha s2 Casein Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Ms
Immunogen Description	KLH conjugated synthetic peptide aa 150-210 222 derived from cow Alpha s2 Casein
Conjugates	HRP
Target Name	Alpha s2 Casein
Other Names	Alpha-S2-casein; CSN1S2
Accession No.	Swiss-Prot#P02663NCBI Gene ID282209
Uniprot	P02663
GeneID	282209;
Excitation Emission	N A
Cell Localization	Secreted
Concentration	1mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Important role in the capacity of milk to transport calcium phosphate. Casocidin-I inhibits the growth of E.coli and S.carnosus.

Note: This product is for in vitro research use only