

Prolactin Rabbit mAb

Catalog No: #56350

Package Size: #56350-1 50ul #56350-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Prolactin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	Prolactin Antibody detects endogenous levels of total Prolactin
Immunogen Description	A synthesized peptide derived from human Prolactin
Other Names	GHA1; Mamotropin; PRL; Prolactin;
Accession No.	Uniprot:P01236
Uniprot	P01236
Calculated MW	23kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000

Images



Western blot analysis of Prolactin expression in Human fetal pituitary lysate.

Product Description

Various hormones are secreted from the anterior pituitary during development and growth. Prolactin is a growth factor secreted by the anterior pituitary that is necessary for the proliferation and differentiation of the mammary glands.

Note: This product is for in vitro research use only