

IGHD Rabbit mAb

Catalog No: #56351

Package Size: #56351-1 50ul #56351-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	IGHD Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	IGHD Antibody detects endogenous levels of total IGH
Immunogen Description	A synthesized peptide derived from human IGH
Other Names	IGHD; IgD; Immunoglobulin delta;
Accession No.	Uniprot:P01880
Uniprot	P01880
Calculated MW	84kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:100~1:500

Images



Western blot analysis of IGH expression in Human tonsil cell lysate.

Product Description

There are two named isoforms of IGH produced by alternative splicing, a secreted form and a single pass type I membrane protein.

Note: This product is for in vitro research use only