

Myogenin Rabbit mAb

Catalog No: #56388

Package Size: #56388-1 50ul #56388-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Myogenin Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC FC
Species Reactivity	Human Mouse Rat
Specificity	Myogenin Antibody detects endogenous levels of total Myogenin
Immunogen Description	A synthesized peptide derived from human Myogenin
Other Names	bHLHc3; cb553; MYF4; MYOG; Myogenic factor 4;
Accession No.	Uniprot:P15173
Uniprot	P15173
Calculated MW	36kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200FC:1:40

Images



Western blot analysis of Myogenin expression in C2C12 cell lysate.

Product Description

Involved in muscle differentiation (myogenic factor). Induces fibroblasts to differentiate into myoblasts. Probable sequence specific DNA-binding protein.

Note: This product is for in vitro research use only