

TAGLN Rabbit mAb

Catalog No: #56409

Package Size: #56409-1 50ul #56409-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TAGLN Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	TAGLN Antibody detects endogenous levels of total TAGLN
Immunogen Description	A synthesized peptide derived from human TAGLN
Other Names	TAGLN; Transgelin; SM22; SM22 alpha; WS3-10;
Accession No.	Uniprot:Q01995
Uniprot	Q01995
Calculated MW	23kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200

Images



Western blot analysis of TAGLN expression in Human stomach lysate.

Product Description

Actin cross-linking/gelling protein (By similarity). Involved in calcium interactions and contractile properties of the cell that may contribute to replicative senescence.

Note: This product is for in vitro research use only