

Parathyroid Hormone Receptor 1 Rabbit mAb

Catalog No: #56417

Package Size: #56417-1 50ul #56417-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Parathyroid Hormone Receptor 1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human
Specificity	Parathyroid Hormone Receptor 1 Antibody detects endogenous levels of total Parathyroid Hormone Receptor 1
Immunogen Description	A synthesized peptide derived from human Parathyroid Hormone Receptor 1
Other Names	PTH receptor; PTH1 receptor; PTH1R; PTHR 1; PTHR; PTHR1;
Accession No.	Uniprot:Q03431
Uniprot	Q03431
Calculated MW	66kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IP:1:50

Images



Western blot analysis of Parathyroid Hormone Receptor 1 expression in HeLa cell lysate.

Product Description

Receptor for parathyroid hormone and for parathyroid hormone-related peptide. The activity of this receptor is mediated by G proteins which activate adenylyl cyclase and also a phosphatidylinositol-calcium second messenger system.

Note: This product is for in vitro research use only