

## CD66b Rabbit mAb

Catalog No: #56442

Package Size: #56442-1 50ul #56442-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	CD66b Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	CD66b Antibody detects endogenous levels of total CD66b
Immunogen Description	A synthesized peptide derived from human CD66b
Other Names	CD66b; CD67; CEACAM8; CGM6; NCA95;
Accession No.	Uniprot:P31997
Uniprot	P31997
Calculated MW	75kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 $\Lambda$ C short term. Store at -20 $\Lambda$ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000

## Images



Western blot analysis of CD66b expression in Human plasma lysate.

## Product Description

Expressed in leukocytes of chronic myeloid Leukemia patients and bone marrow.

---

Note: This product is for in vitro research use only