

RNF40 Rabbit mAb

Catalog No: #56487

Package Size: #56487-1 50ul #56487-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	RNF40 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC IP
Species Reactivity	Mouse, Rat, Human
Specificity	RNF40 Antibody detects endogenous levels of total RNF40
Immunogen Description	A synthesized peptide derived from human RNF40
Other Names	BRE 1B; BRE1 B; BRE1B; RBP95; Rnf40; STARING;
Accession No.	Uniprot:O75150
Uniprot	O75150
Calculated MW	130kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200IP:1:50

Images



Western blot analysis of RNF40 expression in (1) HeLa cell lysate; (2) RAW264.7 cell lysate.

Product Description

E3 ubiquitin-protein ligase that mediates monoubiquitination of 'Lys-120' of histone H2B (H2BK120ub1). H2BK120ub1 gives a specific tag for epigenetic transcriptional activation and is also prerequisite for histone H3 'Lys-4' and 'Lys-79' methylation (H3K4me and H3K79me, respectively).

Note: This product is for in vitro research use only