

NPC1L1 Rabbit mAb

Catalog No: #56507

Package Size: #56507-1 50ul #56507-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NPC1L1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	NPC1L1 Antibody detects endogenous levels of total NPC1L1
Immunogen Description	A synthesized peptide derived from human NPC1L1
Other Names	NPC1 LIKE 1; NPC1L1;
Accession No.	Uniprot:Q9UHC9
Uniprot	Q9UHC9
Calculated MW	149kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000

Images



Western blot analysis of Niemann Pick C1 Like 1 expression in HeLa cell lysate.

Product Description

Play a major role in cholesterol homeostasis. Is critical for the uptake of cholesterol across the plasma membrane of the intestinal enterocyte. Is the direct molecular target of ezetimibe, a drug that inhibits cholesterol absorption.

Note: This product is for in vitro research use only