

GCDFP 15 Rabbit mAb

Catalog No: #56544

Package Size: #56544-1 50ul #56544-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GCDFP 15 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC IP
Species Reactivity	Human
Specificity	GCDFP 15 Antibody detects endogenous levels of total GCDFP 15
Immunogen Description	A synthesized peptide derived from human GCDFP 15
Other Names	BRST2; GCDGP15; gp17; GPIP4; PIP; Prolactin induced protein; SABP; Secretory actin binding protein;
Accession No.	Uniprot:P12273
Uniprot	P12273
Calculated MW	17 kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000ICC/IF:1:50~1:200IP:1:50

Images



Western blot analysis of GCDFP 15 expression in T47D cell lysate.

Product Description

Belongs to the PIP family.

Note: This product is for in vitro research use only