

TNFSF9 Rabbit mAb

Catalog No: #56562

Package Size: #56562-1 50ul #56562-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	TNFSF9 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IP
Species Reactivity	Human
Specificity	TNFSF9 Antibody detects endogenous levels of total TNFSF9
Immunogen Description	A synthesized peptide derived from human TNFSF9
Other Names	4 1BB L; 4 1BBL; 4-1BB ligand; 4-1BBL; Cd137l; Cd157l; Ly63l; Tnfsf9; TNLG5A;
Accession No.	Uniprot:P41273
Uniprot	P41273
Calculated MW	25-32kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IP:1:50

Images



Western blot analysis of TNFSF9 expression in HEK293 cell lysate.

Product Description

Cytokine that binds to TNFRSF9. Induces the proliferation of activated peripheral blood T-cells. May have a role in activation-induced cell death (AICD). May play a role in cognate interactions between T-cells and B-cells/macrophages.

Note: This product is for in vitro research use only