

TEAD1/2/3/4 Rabbit mAb

Catalog No: #56565

Package Size: #56565-1 50ul #56565-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TEAD1/2/3/4 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB
Species Reactivity	Human
Specificity	TEAD1 + TEAD2 + TEAD3 + TEAD4 Antibody detects endogenous levels of total TEAD1 + TEAD2 + TEAD3 + TEAD4
Immunogen Description	A synthesized peptide derived from human TEAD1 + TEAD2 + TEAD3 + TEAD4
Other Names	DTEF 1; NTEF 1; TCF13; TEAD1; TEAD2; TEAD3; TEAD4; TEAD5; TEF1; TEF4; TEF5;
Accession No.	Uniprot:P28347/Q15561/Q15562/Q99594
Uniprot	P28347/Q15561/Q15562/Q99594
Calculated MW	54kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:1000

Images



Western blot analysis of TEAD1/2/3/4 expression in SW480 cell lysate.

Product Description

Transcription factor which plays a key role in the Hippo signaling pathway, a pathway involved in organ size control and tumor suppression by restricting proliferation and promoting apoptosis.

Note: This product is for in vitro research use only