

NUDC Rabbit mAb

Catalog No: #56570

Package Size: #56570-1 50ul #56570-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	NUDC Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP FC
Species Reactivity	Human Mouse Rat
Specificity	NUDC Antibody detects endogenous levels of total NUDC
Immunogen Description	A synthesized peptide derived from human NUDC
Other Names	HNUDC; MNUDC; NPD011; nudC; SIG92;
Accession No.	Uniprot:Q9Y266
Uniprot	Q9Y266
Calculated MW	42kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200ICC/IF:1:50~1:200IP:1:50FC:1:50

Images



Western blot analysis of NUDC expression in HeLa cell lysate.

Product Description

Plays a role in neurogenesis and neuronal migration (By similarity). Necessary for correct formation of mitotic spindles and chromosome separation during mitosis. Necessary for cytokinesis and cell proliferation.

Note: This product is for in vitro research use only