

NAP1L1 Rabbit mAb

Catalog No: #56643

Package Size: #56643-1 50ul #56643-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	NAP1L1 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB FC
Species Reactivity	Human Mouse Rat
Specificity	NAP1L1 Antibody detects endogenous levels of total NAP1L1
Immunogen Description	A synthesized peptide derived from human NAP1L1
Other Names	NRP; NAP1; NAP1L;
Accession No.	Uniprot:P55209
Uniprot	P55209
Calculated MW	52,60kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000FC:1:50

Images



Western blot analysis of NAP1L1 expression in Jurkat cell lysate.

Product Description

May be involved in modulating chromatin formation and contribute to regulation of cell proliferation.

Note: This product is for in vitro research use only