

CIRBP Rabbit mAb

Catalog No: #56906



Package Size: #56906-1 50ul #56906-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	CIRBP Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC ICC/IF IP
Species Reactivity	Human
Specificity	CIRBP Antibody detects endogenous levels of total CIRBP
Immunogen Description	A synthesized peptide derived from CIRBP
Other Names	A18HNRNP; cirbp; CIRP; Cold inducible RNA binding protein;
Accession No.	Uniprot:Q14011
Uniprot	Q14011
Calculated MW	19kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4Λ C short term. Store at -20Λ C long term. Avoid freeze / thaw cycle.


Application Details

WB:1:500~1:2000
IHC:1:50~1:200
ICC/IF:1:50~1:200
IP:1:50

Images



Western blot analysis of CIRBP expression in HepG2 cell lysate.



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.

Product Description

Cold-inducible mRNA binding protein that plays a protective role in the genotoxic stress response by stabilizing transcripts of genes involved in cell survival. Acts as a translational activator.

Note: This product is for in vitro research use only