

CD84 Rabbit mAb

Catalog No: #56987

Package Size: #56987-1 50ul #56987-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CD84 Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB IHC
Species Reactivity	Human
Specificity	CD84 Antibody detects endogenous levels of total CD84
Immunogen Description	A synthesized peptide derived from human CD84
Other Names	CD84; CDw84; hCD84; LY9B; mCD84; SLAMF5;
Accession No.	Uniprot:Q9UIB8
Uniprot	Q9UIB8
Calculated MW	60-85kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4 Λ C short term. Store at -20 Λ C long term. Avoid freeze / thaw cycle.

Application Details

WB:1:500~1:2000IHC:1:50~1:200

Images



Western blot analysis of CD84 expression in Raji cell lysate.

Product Description

Plays a role as adhesion receptor functioning by homophilic interactions and by clustering. Recruits SH2 domain-containing proteins SH2D1A/SAP.

Note: This product is for in vitro research use only