

## PPP6C Rabbit mAb

Catalog No: #56998

Package Size: #56998-1 50ul #56998-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	PPP6C Rabbit mAb
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB ICC/IF FC
Species Reactivity	Human
Specificity	PPP6C Antibody detects endogenous levels of total PPP6C
Immunogen Description	A synthesized peptide derived from human PPP6C
Other Names	PP6; PP6C; PPP6; Ppp6c;
Accession No.	Uniprot:O00743
Uniprot	O00743
Calculated MW	32kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4ℳ C short term. Store at -20ℳ C long term. Avoid freeze / thaw cycle.

## Application Details

WB:1:500~1:2000ICC/IF:1:50~1:200FC:1:50

## Images



Western blot analysis of PPP6C expression in HeLa cell lysate.

## Product Description

Catalytic subunit of protein phosphatase 6 (PP6). PP6 is a component of a signaling pathway regulating cell cycle progression in response to IL2 receptor stimulation.

---

Note: This product is for in vitro research use only