

SP-C Antibody PE Conjugated

Catalog No: #C00573P

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SP-C Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 24-58 197 derived from human SP-C
Conjugates	PE
Target Name	SP-C
Other Names	SP-C; PSP-C; SFTP2; SMDP2; BRICD6; Pulmonary surfactant-associated protein C; Pulmonary surfactant-associated proteolipid SPL(Val); SP5; SFTPC
Accession No.	Swiss-Prot#P11686NCBI Gene ID6440
Uniprot	P11686
GeneID	6440;
Excitation Emission	480,565nm 578nm
Cell Localization	Extracellular, Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Background

Pulmonary surfactant associated proteins promote alveolar stability by lowering the surface tension at the air-liquid interface in the peripheral air spaces.

Note: This product is for in vitro research use only