

Tyrosine Hydrolase Mouse mAb

Catalog No: #44222

Package Size: #44222-1 50ul #44222-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Tyrosine Hydrolase Mouse mAb
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Applications	IHC-p;WB
Species Reactivity	Human; Mouse;Rat
Specificity	This antibody detects endogenous levels of Tyrosine Hydrolase at Human, Mouse,Rat
Immunogen Description	Synthesized peptide derived from human Tyrosine Hydrolase
Accession No.	Swiss Prot:P07101GeneID:7054
Calculated MW	55kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.30% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:1000-2000

IHC-p1:50-200

Images



Whole cell lysates of MDCK were separated by 10% SDS-PAGE, and the membrane was blotted with anti-Tyrosine Hydroxylase antibody.



Immunohistochemical analysis of paraffin-embedded human Colon cancer.

Note: This product is for in vitro research use only