

# Anti-CD107a Antibody

Catalog No: #44226



Package Size: #44226-1 50ul #44226-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Anti-CD107a Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	This antibody is purified through a protein G column.
Applications	WB, IHC, IF/ICC
Species Reactivity	Human, Mk
Immunogen Description	Recombinant fusion protein of human CD107a.
Other Names	Lysosome-associated membrane glycoprotein 1; LAMP-1; Lysosome-associated membrane protein 1; CD107 antigen-like family member A; CD107a
Accession No.	Swiss Prot:P11279GenelD:3916
Formulation	Mouse IgG1 kappa. Liquid in PBS, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

WB (1/500 - 1/1000),

IHC (1/50 - 1/200), IF/ICC (1/10 - 1/50)

Note: This product is for in vitro research use only